



α -Smooth muscle actin (α -SM1; 1A4)

CATALOG NUMBER: MUB0100S

CLONE: 1A4

SPECIES / ISOTYPE: mouse IgG2a

PRODUCT FORM: culture supernatant of monoclonal antibody

BACKGROUND

Among the six actin isoforms described in mammals, two are found in virtually all cells (β - and γ -cytoplasmic), two are detected in smooth muscle cells (α - and γ -smooth muscle) and two are present in striated muscles, one predominantly in skeletal (α -skeletal) and one in cardiac (α -cardiac) muscle cells. These actin isoforms differ slightly in their N-terminus and the sequences of each of them are perfectly conserved in higher vertebrates. Alpha-smooth muscle actin is abundant in vascular and visceral smooth muscle cells. In addition, it has also been shown to appear in stress fibers of fibroblastic cells during pathological situations involving contractile phenomena such as wound healing and fibrocontractive diseases.

SOURCE

α -SM1 (clone 1A4) is a mouse monoclonal IgG2a antibody derived by fusion of Sp 2/0 mouse myeloma cells with spleen cells from a BALB/c mouse immunized with a peptide comprising the first 10 amino acids of α -smooth muscle actin with an acetylated N-terminus coupled to keyhole limpet hemocyanin via the C-terminal cysteine (Ac-EEEDSTALVC).

PRODUCT

Each vial contains 1ml of culture supernatant monoclonal antibody containing 0.09% sodium azide.

SPECIFICITY

α -SM1 reacts exclusively with α -smooth muscle actin which is typical for vascular and visceral smooth muscle cells, but which is also present in myofibroblasts. The epitope that is recognized by α -SM1 is Ac-EEED.

α -SM1 is useful for immunocytochemistry, immunohistochemistry on frozen and paraffin-embedded tissues, immunoblotting, immunoelectron microscopy and ELISA. Optimal antibody dilution should be determined by titration; recommended range is 1:10 – 1:100 for immunohistochemistry with avidin-biotinylated

horseradish peroxidase complex (ABC) as detection reagent, and 1:100 – 1:500 for immunoblotting applications.

SPECIES REACTIVITY

The epitope recognized by α -SM1 is highly conserved. The antibody therefore cross-reacts with many species including protochordates, lower craniates and mammals (amongst others human, rat, chicken, goat, rabbit, mouse, swine).

STORAGE

Store at 4°C, **do not freeze.**

REFERENCES

1. Skalli, O., Ropraz, P., Trzeciak, A., Benzonana, G., Gillesse, D., and Gabbiani, G. (1986). A monoclonal antibody against alpha-smooth muscle actin: a new probe for smooth muscle differentiation. *J Cell Biol* 103, 2787-96.
2. Skalli, O., Schurch, W., Seemayer, T., Lagace, R., Montandon, D., Pittet, B., and Gabbiani, G. (1989). Myofibroblasts from diverse pathologic settings are heterogeneous in their content of actin isoforms and intermediate filament proteins. *Lab Invest* 60, 275-85.
3. Babaev, V.R., Bobryshev, Y.V., Stenina, O.V., Tararak, E.M. and Gabani, G. (1990). Heterogeneity of smooth muscle in atherosclerotic plaque of human aorta. *American Journal of Pathology* 136, 1031-42.
4. Sappino, A. P., Schurch, W., and Gabbiani, G. (1990). Differentiation repertoire of fibroblastic cells: expression of cytoskeletal proteins as marker of phenotypic modulations. *Lab Invest* 63, 144-61.
5. Vyalov, S. L., Gabbiani, G., and Kapanci, Y. (1993). Rat alveolar myofibroblasts acquire alpha-smooth muscle actin expression during bleomycin-induced pulmonary fibrosis. *Am J Pathol* 143, 1754-65.
6. Chaponnier, C., Goethals, M., Janmey, P. A., Gabbiani, F., Gabbiani, G., and Vandekerckhove, J. (1995). The specific NH2-terminal sequence Ac-EEED of alpha-smooth

WARNING and CAUTION

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

This product contains sodium azide. To prevent formation of toxic vapors, do not mix with strong acidic solutions. To prevent formation of potentially explosive metallic azides in metal plumbing, always wash into drain with copious quantities of water.

This datasheet is as accurate as reasonably achievable, but MUBio Products BV accepts no liability for any inaccuracies or omissions in this information.



- muscle actin plays a role in polymerization in vitro and in vivo. *J Cell Biol* **130**, 887-95.
7. Simoncelli, F., Fagotti, A., Di Rosa, I., Panara, F., Chaponnier, C., Gabbiani, G., Pascolini, R., (1996). Expression of an actin in protochordates and lower craniates defined by anti- α SM-1. *European J of Cell Biol* **69**, 297-300.
 8. Hinz, B., Celetta, G., Tomasek, J. J., Gabbiani, G., and Chaponnier, C. (2001). Alpha-smooth muscle actin expression upregulates fibroblast contractile activity. *Mol Biol Cell* **12**, 2730-41.
 9. Hinz, B., Gabbiani, G., and Chaponnier, C. (2002). The NH2-terminal peptide of alpha-smooth muscle actin inhibits force generation by the myofibroblast in vitro and in vivo. *J Cell Biol* **157**, 657-63.
 10. Hinz, B., Dugina, V., Ballestrem, C., Wehrle-Haller, B., Chaponnier, C. (2003). Alpha-smooth muscle actin is crucial for focal adhesion maturation in myofibroblasts. *Mol Biol Cell* **14**, 2508-19.
 11. Clément, S., Hinz, B., Dugina, V., Gabbiani, G. and Chaponnier, C. (2004). The N-terminal Ac-EEED sequence plays a role in alpha-smooth-muscle actin incorporation into stress fibers. *Journal of Cell Science* **118**, 1395-1404.
 12. Chaponnier, C. and Gabbiani, G. (2004). Pathological situation characterized by altered actin isoform expression. *J Pathol* **204**, 386-95.
 13. Clément, S., Stouffs, M., Bettiol, E., Kampf, S., Krause, K., Chaponnier, C. and Jaconi, M. (2006). Expression and function of alpha-smooth muscle actin during embryonic-stem-cell-derived cardiomyocyte differentiation. *Journal of Cell Science* **120**, 229-38.

© 2010 MUBio Products B.V.
Datasheet version: MUB_0100P_100202

WARNING and CAUTION

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

This product contains sodium azide. To prevent formation of toxic vapors, do not mix with strong acidic solutions. To prevent formation of potentially explosive metallic azides in metal plumbing, always wash into drain with copious quantities of water.

This datasheet is as accurate as reasonably achievable, but MUBio Products BV accepts no liability for any inaccuracies or omissions in this information.